D.B.M. C-LOCK 4%

STERILE SODIUM CITRATE 4%

Sodium citrate 4% sterile "lock" solution for intravenous catheters.
 Sterile - Latex-free - Disposable

QUALITATIVE AND QUANTITATIVE COMPOSITION
 SODIUM CITRATE TRIBASIC DIHYDRATE 40 g
 Excipients: ACIDO CITRICO MONOHYDRATE q.s. to pH 6.7 - WFI q.s. to 1000 mL

DEVICE TYPE
 Sterile, pyrogen-free, ready-to-use Medical Device - Class Ib.
 Medical Device not subject to a medical prescription

- PACKAGING**
PC0131A1: 3 mL disposable pre-filled syringe, filled to 2.5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A2: 3 mL disposable pre-filled syringe, filled to 2.5 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
PC0131A3: 3 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A4: 3 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A5: 3 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
PC0131A6: 5 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A7: 5 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
PC0131A8: 5 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A9: 10 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A10: 10 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A11: 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A12: 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch and in a second lumen pouch.
PC0131A13: 20 mL disposable pre-filled syringe, filled to 15 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PC0131A14: 20 mL disposable pre-filled syringe, filled to 20 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.

Only the contents and the outside of the syringe are sterile with all types of packaging. The syringe and the product are moist-heat stable in autoclave.

The volume reported on the label refers to the filling volume of the syringe. The graduated scale reported on the label of the empty syringe (without syringe) has no measuring function.

INTENDED USE
 The Medical Device is intended for the safe practice of vascular access. Its function is to maintain the central venous catheters (CVC) patent in the intertidal interval, also performing an anticoagulant function. The presence of citrate triggers a chelating action for the calcium ions in the blood and tissues, making them no longer available for the coagulative cascade. Therefore, this solution eliminates the risks of systemic anticoagulation connected with the use of heparin as a lock solution and ensures total safety and handiness thanks to the ready-to-use pre-filled syringe.

By adopting an aseptic procedure, D.B.M. C-LOCK 4% can be used on the sterile field.

SPECIAL INSTRUCTIONS AND PRECAUTIONARY MEASURES

- To be used only by suitably trained medical staff who are capable to take all necessary precautions.
- Do not administer intravenously directly.
- Do not use the product for purposes other than its intended use.
- Do not use the product if the packaging is open or damaged. If the syringe is not perfectly intact or if the cap is so damaged that the syringe leaks.
- Do not use the product if the solution is not colorless or if it contains particulate matter.
- Once opened, the device must be used immediately and disposed of after use.
- Single-use only. Dispose any partially used solution in a way which may cause cross contamination between patient and operator.
- Do not re-sterilize.
- Do not let air into the syringe and do not allow air to get trapped in the fluid path. Use only with compatible Luer Lock adapters.
- If the vascular access device is not adequately sealed, closing may result in thrombosis.
- Do not use an injection volume greater than the lumen primary volume specified by the manufacturer of the vascular access device.
- Should blood be taken through the catheter, consider that even minimal residues of sodium citrate can lead to incorrect clinical results.
- Keep out of reach and sight of children.

CONTRAINDICATIONS AND SIDE EFFECTS

If a portion of sodium citrate, which is inside the catheter during the intertidal interval, inadvertently gets into the patient's bloodstream, chelation of calcium ions may occur, resulting in hypocalcemia and hypernatremia, as well as possible emergence of parasthesia and/or dysesthesia.


INSTRUCTIONS FOR USE

- Check the validity of the expiry date reported on the sealing side of the pouch.
- Open the pouch and take out the syringe. Check the correct position of the syringe cap, the transparency of the solution and the absence of particulate matter. If necessary, place the syringe on a sterile field.
- Clear the catheter from any previously used closure solutions by sucking up with an empty sterile syringe.
- Wash the lumens of each chamber, using a sterile isotonic solution (Sodium Chloride 0.9%), to clean up any residual haematic components from previous operations.
- Flush the chamber, with closure cap inserted, to remove the rubber seal.
- Remove the protective cap, making sure that no contact contamination of the syringe Luer cone occurs.
- Remove the syringe in a vertical position and eliminate any air bubbles.
- Connect the syringe to the Luer Lock adapter on the vascular access device while taking care to avoid any contact contamination of the syringe Luer cone.
- Slowly inject (5-10 seconds) the amount of "lock" solution in accordance with the instructions of the manufacturer of the vascular access device and according to the procedure adopted by your institution.
- Do not inject an amount exceeding the volume of the catheter lumen.
- After use, dispose of the empty syringe and any unused solution volume in accordance with the national regulations in force.

SHELF-LIFE
 3 years in an intact pouch.

STORAGE

Store the syringes tightly sealed, away from direct light and heat sources, at a temperature between +5 °C and +30 °C. Do not freeze. The expiration date refers to the product in its intact packaging and correctly stored. Do not use the Medical Device after the expiration date.

SYMBOLS USED ON THE LABEL	
	Do not reuse - Disposable
	Manufacturer
	Product code
	Moist-heat sterilized
	Production batch
	Expiry date
	Temperature limits for storage
	Refer to the instructions for use
	Latex-free
 Medical Device compliant with Directive 93/42/CEE	
D.B.M. Srl Via dell'Industria, 2 - 23030 Tovo di Sant'Agata (SO) - Italy T +39 0342 771613 - info@dbmitalia.it - www.dbmitalia.it	

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