	0476 D.B.M. S-GEL	
STERILE LUBRICATING GEL		
Pre-filled syringe with sterile lubricating gel Sterile – Latex-free – Disposable		

QUALITATIVE AND QUANTITATIVE COMPOSITION

GELLED HYDROPOLYHYDROSILIC SOLUTION.

Excipients: demineralized water, carboxyvinyl polymer, sodium hydroxide, methyl P-oxbenzoate, propylene glycol, hydroxypropyl methylcellulose.

DEVICE TYPE

Sterile, ready-to-use medical device - Class IIa.
Medical Device not subject to a medical prescription.

PACKAGING

PCF09MA1: 5 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PCG09MA1: 10 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PCG09MA1: 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
PAF09MA1: 5 mL disposable pre-filled syringe, filled to 5 mL, with Luer Slip adapter, individually packaged in a PET/PP paper pouch.
PAF09MA1: 10 mL disposable pre-filled syringe, filled to 5 mL, with Luer Slip adapter, individually packaged in a PET/PP paper pouch.
PAF09MA1: 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Slip adapter, individually packaged in a PET/PP paper pouch.
PAF09NE1: 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Slip adapter, individually packaged in a PET/PP paper pouch and in a second aluminium pouch.

Only the contents and the outside of the syringe are sterile with all types of packaging. The syringe and the product are moist-heat sterilized in autoclave.

The volume reported on the label refers to the filling volume of the syringe. The graduated scale reported on the label of the primary packaging (syringe) has no measuring function.

INTENDED USE

This Medical Device is intended to be used as a lubricant on endoscopes, Foley catheters, specula, endotracheal tubes, endotracheal tubes, etc. used in diagnostic and therapeutic procedures, in order to facilitate their insertion by reducing friction between the contact surfaces (tissues) and facilitating their sliding. It can also be used in ultrasound inspections (ultrasonography) as a coupling medium and for the preparation of mixtures for diagnostic and therapeutic purposes, e.g. mixtures containing contrast medium for retrograde radiological visualization of the urethral or ureteral profile.

By adopting an aseptic procedure, D.B.M. S-Gel can be used on the sterile field.

SPECIAL INSTRUCTIONS AND PRECAUTIONARY MEASURES

- Do not administer intravenously directly. Do not use the product for purposes other than its intended use.
- The gel is not electrically conductive.
- Do not apply the gel to skin with wounds or erythema, or to patients with mucous membrane damage or bleeding.
- Do not use the product if the packaging is open or damaged, if the syringe is not perfectly intact or if the cap is so damaged that the syringe leaks.
- Do not use the product if the solution is not colourless or if it contains particulate matter.
- Once opened, the device must be used immediately and disposed of after use.
- Single-use only. Disposal of any partially used product, the reuse of which may cause cross contamination between patient and operator.
- Do not re-sterilize.
- Keep out of reach and sight of children.

CONTRAINDICATIONS AND SIDE EFFECTS

This Medical Device may cause irritation to allergic patients. Do not use in patients with hypersensitivity to alkyl-4-hydroxybenzoates (parabens) or other components. In case of hypersensitivity, skin reactions may show up, with redness and burning.

INSTRUCTIONS FOR USE

- In order to ensure a sterile application of the product, the health-care worker must comply with the procedures provided for by the health-care institution, applying an aseptic technique to the entire procedure.
- Check the validity of the expiry date reported on the sealing side of the pouch.
- Open the pouch and take out the syringe. Check the correct position of the syringe cap and the transparency of the solution. If necessary, place the syringe on the sterile field.
- Remove the protective cap, making sure that no contact contamination of the syringe Luer cone occurs.
- Apply to the skin, mucosa or instrument to be lubricated (e.g. catheter, tube or endoscope), according to the procedure adopted by your health-care institution.
- After use, dispose of the empty syringe and any unused gel volumes in accordance with the national regulations in force.

SHELF-LIFE

3 years in an intact pouch.

STORAGE

Store the syringes tightly sealed, away from direct light and heat sources, at a temperature between + 5 °C and + 30 °C. Do not freeze. The expiration date refers to the product in its intact packaging and correctly stored. Do not use the Medical Device after the expiration date.

SYMBOLS USED ON THE LABEL	
	Do not reuse - Disposable
	Manufacturer
	Product code
	Moist-heat sterilized
	Production batch
	The outside of the secondary packaging (PET/PP paper pouch) is not sterile
	Do not re-sterilize.
	Temperature limits for storage
	Keep away from light and heat sources
	Latex-free
	Protect from humidity
	Medical Device compliant with Directive 93/42/CEE
	D.B.M. Srl Via dell'industria, 2 - 23030 Tovo di Sant'Agata (SO) - Italy T +39 0342 771613 - info@dbmitalia.it - www.dbmitalia.it

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