<b>CE</b> 0476 <b>D.B.M. C-LOCK 30%</b>	<b>EN</b>
STERILE SODIUM CITRATE 30%	
Sodium citrate 30% sterile "lock" solution for intravenous catheters. Sterile - Latex-free - Disposable	
QUALITATIVE AND QUANTITATIVE COMPOSITION SODIUM CITRATE TRIBASIC DIHYDRATE 300 g Expipients: CITRIC ACID MONOHYDRATE q.s. pH 6.7-5 - WFI q.s. to 1000 mL	
DEVICE TYPE Sterile, pyrogen-free, ready-to-use Medical Device - Class IIB. Medical Device not subject to a medical prescription	

QUALITATIVE AND QUANTITATIVE COMPOSITION  
SODIUM CITRATE TRIBASIC DIHYDRATE 300 g  
Expipients: CITRIC ACID MONOHYDRATE q.s. pH 6.7-5 - WFI q.s. to 1000 mL

DEVICE TYPE  
Sterile, pyrogen-free, ready-to-use Medical Device - Class IIB.  
Medical Device not subject to a medical prescription

- PACKAGING**
- PC15K1AZ:** 3 ml disposable pre-filled syringe, filled to 2.5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15K2AZ:** 3 ml disposable pre-filled syringe, filled to 2.5 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15L1AZ:** 3 ml disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15L2AZ:** 3 ml disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15L1LAZ:** 3 ml disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15L2LAZ:** 3 ml disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15L1LAZ:** 5 ml disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15L2LAZ:** 5 ml disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15L1LAZ:** 10 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15L2LAZ:** 10 mL disposable pre-filled syringe, filled to 3 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15M1AZ:** 10 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15M2AZ:** 10 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15L1LAZ:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15L2LAZ:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.
  - PC15M1AZ:** 20 mL disposable pre-filled syringe, filled to 15 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.
  - PC15M2AZ:** 20 mL disposable pre-filled syringe, filled to 15 mL, with Luer Lock adapter, twin packaged in a PET/PP paper pouch.

Only the contents and the outside of the syringe are sterile with all types of packaging. The syringe and the product are most-heat sterilized in autoclave.

The volume reported on the label refers to the filling volume of the syringe. The graduated scale reported on the label of the primary packaging (syringe) has no measuring function.

**INTENDED USE**

This Medical Device is intended for the safe practice of vascular access. Its function is to maintain the central venous catheter (CVC) patent in the intertidal interval, also performing an anticoagulant function. The presence of citrate triggers a chelating action for the calcium ions in the blood and tissues, making them no longer available for the coagulative cascade. Therefore, this solution eliminates the risks of systemic anticoagulation connected with the use of heparin as a lock solution and ensures total safety and hardness thanks to the ready-to-use pre-filled syringe.

D.B.M. C-LOCK 30% also acts as an antimicrobial effect and reduces biofilm formation in vascular access devices by binding and removing calcium and magnesium ions from the surrounding environment.

By adopting an aseptic procedure, D.B.M. C-LOCK 30% can be used on the sterile field.

**SPECIAL INSTRUCTIONS AND PRECAUTIONARY MEASURES**

- To be used only by suitably trained medical staff who are capable to take all necessary precautions.
- Do not administer intravenously directly.
- Do not use the product for purposes other than its intended use.
- Do not use the product if the packaging is open or damaged, if the syringe is not perfectly intact or if the cap is so damaged that the syringe leaks.
- Do not use the product if the solution is not colorless or if it contains particulate matter.
- Once opened, the device must be used immediately and disposed of after use.
- Single-use only. Dispose of any partially used product, the reuse of which may cause cross contamination between patient and operator.
- Do not re-sterilize.
- Do not immerse the syringe and do not allow air to get trapped in the fluid path. Use only with compatible Luer Lock adapters.
- The vascular access device is not adequately sealed, oils may form, resulting in thrombosis.
- Do not use an injection volume greater than the lumen priming volume specified by the manufacturer of the vascular access device.
- Should blood be taken through the catheter, consider that even minimal residues of sodium citrate can lead to incorrect clinical results.
- Keep out of reach and sight of children.

**CONTRAINDICATIONS AND SIDE EFFECTS**

If a portion of sodium citrate, which is inside the catheter during the intertidal interval, inadvertently gets into the patient's bloodstream, chelation of calcium ions may occur, resulting in hypocalcemia and hypervolemia, as well as possible emergence of parathesias and/or dysesthesia.

**INSTRUCTIONS FOR USE**

- Disposable device for a single patient.
- In order to ensure a sterile application of the product, the health-care worker must comply with the procedures provided for by the health-care institution, applying an aseptic technique to the entire procedure.
- Check the validity of the expiry date reported on the sealing side of the pouch.
- Open the pouch and take out the syringe. Check the correct position of the syringe cap, the transparency of the solution and the absence of particulate matter. If necessary, place the syringe on a sterile field.
- Clear the catheter from any previously used closure solutions by sucking up with an empty sterile syringe.
- Wash the lumens of each catheter, using a sterile isotonic solution (Sodium Chloride 0.9%), to clean up any residual haematic components from previous operations.
- With the closure cap inserted, push the plunger to remove the rubber seal.
- Remove the protective cap, making sure that no contact contamination of the syringe Luer cone occurs.
- Keep the syringe in a vertical position and eliminate any air bubbles.
- Connect the syringe to the Luer Lock adapter. On the vascular access device avoid taking care to avoid any contact contamination of the syringe Luer cone.
- Slowly inject (5-10 seconds) the amount of "lock" solution in accordance with the instructions of the manufacturer of the vascular access device and according to the procedure adopted by your institution.
- Do not inject an amount exceeding the volume of the catheter specified by the manufacturer of the catheter.
- After use, dispose of the empty syringe and any unused solution volumes in accordance with the national regulations in force.

**SHELF-LIFE**

3 years in an intact pouch.

**STORAGE**

Store the syringes tightly sealed, away from direct light and heat sources, at a temperature between + 5 °C and + 30 °C. Do not freeze. The expiration date refers to the product in its intact packaging and correctly stored. Do not use the Medical Device after the expiration date.

SYMBOLS USED ON THE LABEL	
	Do not reuse - Disposable
	Product code
	Production batch
	Expiry date
	Do not use if the pouch is open or damaged
	Refer to the instructions for use
	Latex-free
	Medical Device compliant with Directive 93/42/CEE
	Via dell'Industria, 2 - 23030 Tovo di Sant'Agata (SO) - Italy T +39 0342 771613 - info@dbmitalia.it - www.dbmitalia.it

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