**D.B.M. AQUAFLUSH****[EN]****STERILE WATER**

Pre-filled syringe with sterile water for Foley catheters  
Sterile – Latex-free – Disposable

**QUALITATIVE AND QUANTITATIVE COMPOSITION**  
NON-INJECTABLE STERILE WATER.**DEVICE TYPE**

Sterile, ready-to-use medical device - Class IIa.  
Medical Device not subject to a medical prescription

**PACKAGING**  
**PCF01MA1:** 5 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.  
**PCG01MA1:** 10 mL disposable pre-filled syringe, filled to 5 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.  
**PCG01NA1:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch.  
**PCG01NE1:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Lock adapter, individually packaged in a PET/PP paper pouch and in a second aluminium pouch.  
**PAF01MA1:** 5 mL disposable pre-filled syringe, filled to 5 mL, with Luer Slip fitting, individually packaged in a PET/PP paper pouch.  
**PAG01MA1:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Slip fitting, individually packaged in a PET/PP paper pouch.  
**PAG01NA1:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Slip fitting, individually packaged in a PET/PP paper pouch.  
**PAG01NE1:** 10 mL disposable pre-filled syringe, filled to 10 mL, with Luer Slip adapter, individually packaged in a PET/PP paper pouch and in a second aluminium pouch.

Only the contents and the outside of the syringe are sterile with all types of packaging. The syringe and the product are moist-heat sterilized in autoclave.

The volume reported on the label refers to the filling volume of the syringe. The graduated scale reported on the label of the primary packaging (syringe) has no measuring function.

**INTENDED USE**

This Medical Device is exclusively intended to be used for balloon inflation in Foley catheters.

By adopting an aseptic procedure, D.B.M. Aquaflush can be used on the sterile field.

**SPECIAL INSTRUCTIONS AND PRECAUTIONARY MEASURES**

- Do not administer intravenously directly. Do not use the product for purposes other than its intended use.
- Do not use the product if the packaging is open or damaged, if the syringe is not perfectly intact or if the cap is so damaged that the syringe leaks.
- Do not use the product if the solution is not colourless or if it contains particulate matter.
- Once opened, the device must be used immediately and disposed of after use.
- Single-use only. Dispose of any partially used product, the reuse of which may cause cross contamination between patient and operator.
- Do not re-sterilize.
- Keep out of reach and sight of children.

**CONTRAINDICATIONS AND SIDE EFFECTS**

No contraindications. No side effects are known if the solution is used according to its intended use.

**INSTRUCTIONS FOR USE**

Disposable device for a single patient.

- In order to ensure a sterile application of the product, the health-care worker must comply with the procedures provided for by the health-care institution, applying an aseptic technique to the entire procedure.
- Check the validity of the expiry date reported on the sealing side of the pouch.
  - Open the pouch and take out the syringe. Check the correct position of the syringe cap and the transparency of the solution. If necessary, place the syringe on the sterile field.
  - Remove the protective cap, making sure that no contact contamination of the syringe Luer cone occurs.
  - Connect the syringe to the balloon inflation channel on the catheter, taking care to avoid any contamination through contact with the Luer cone of the syringe.
  - Slowly inject the amount of solution for a volume suitable to the correct and effective anchoring of the balloon in the Foley catheter, according to the procedure adopted by your institution.
  - After use, dispose of the empty syringe and any unused solution volumes in accordance with the national regulations in force.

**SHELF-LIFE**

3 years in an intact pouch.

**STORAGE**

Store the syringes tightly sealed, away from direct light and heat sources, at a temperature between + 5 °C and + 30 °C. Do not freeze. The expiration date refers to the product in its intact packaging and correctly stored. Do not use the Medical Device after the expiration date.

SYMBOLS USED ON THE LABEL	
	Do not reuse - Disposable
	Manufacturer
	Product code
	Moist-heat sterilized
	Production batch
	The outside of the secondary packaging (PET/PP paper pouch) is not sterile
	Expiry date
	Do not re-sterilize.
	Temperature limits for storage
	Refer to the instructions for use
	Keep away from light and heat sources
	Latex-free
	Protect from humidity
	Medical Device compliant with Directive 93/42/CEE
	Via dell'Industria, 2 - 23030 Tovo di Sant'Agata (SO) - Italy T +39 0342 771613 - info@dbmitalia.it - www.dbmitalia.it

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