D.B.M. AQUAFUSH

EN
 Syringe prefilled with sterile water for bladder catheters
 Sterile - Latex-free - For single use

STERILE WATER
 Read these instructions for use carefully before each use and keep them easily accessible for interested users.
 Please refer to the specific information on intended use, indication and contraindications.
 Please read these instructions for use carefully. Improper use of the products can cause serious injury to the patient. The user or third parties.
 D.B.M. S.r.l. disclaims all liability towards private persons or persons in the event that manoeuvres other than those described in this instruction sheet are carried out and in cases where the device is not used by qualified personnel.

QUALITATIVE AND QUANTITATIVE COMPOSITION
 Ph. Eur.-compliant WATER

TYPE OF DEVICE
 Sterile, ready-to-use medical device.

AVAILABLE VARIANTS

STERILIZED USING MOIST HEAT	STERILIZED BY GAMMA RADIATION
5 mL pre-filled single-use syringe, filled to 5 mL, with Luer Lock fitting, individually packaged in a paper-PET/PP pouch.	5 mL pre-filled single-use syringe, filled to 5 mL, with Luer Lock fitting, individually packaged in a paper-PET/PE pouch.
5 mL pre-filled single-use syringe, filled to 5 mL, with Luer Lock fitting, individually packaged in a paper-PET/PP pouch and in a second aluminium pouch.	5 mL pre-filled single-use syringe, filled to 5 mL, with Luer Lock fitting, individually packaged in an aluminium-PET/PE pouch.
10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in a paper-PET/PP pouch.	10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in a paper-PET/PE pouch.
10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in a paper-PET/PP pouch and in a second aluminium pouch.	10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in an aluminium-PET/PE pouch.
10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in a paper-PET/PP pouch.	10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in a paper-PET/PE pouch.
10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in a paper-PET/PP pouch and in a second aluminium pouch.	10 mL pre-filled single-use syringe, filled to 10 mL, with Luer Lock fitting, individually packaged in an aluminium-PET/PE pouch.

The volume indicated on the label refers to the syringe filling volume. The graduated scale on the label of the immediate packaging (syringe) has no measuring function.

INTENDED USE AND CLINICAL BENEFIT
 This medical device is intended to be used exclusively for infusing bladder catheter balloons.
 Aqueous solutions have the advantage, unlike saline solutions, that they do not crystallize, making the balloon porous and causing it to deflate with the consequent risk of infection.
 This allows the balloon fill volume to remain constant, ensuring greater safety for catheterized patients. The fact that it is not necessary to check and restore balloon inflation periodically saves healthcare facilities time and money compared to water-only solutions.
 Adopting an aseptic technique, D.B.M. Aquafush can be used in sterile fields.
 It can be used on adult and paediatric patients.

INDICATIONS FOR USE
 This medical device is indicated exclusively for infusing bladder catheter balloons.

PRECAUTIONS
 Do not use the product for purposes other than its intended use.
 Do not use the product if the packaging is open or damaged, if the syringe is not fully intact or if the cap is damaged in such a way as to cause leakage from the syringe.
 Do not use the product if the solution is not colourless or contains visible particulate matter.
 Do not open the device until it is ready to be used.
 For single use only. Dispense a portion of any partially used product, the reuse of which could cause cross contamination between patient and operator.
 Do not draw air into the syringe and do not allow air to be trapped in the fluid path.
 Use only with compatible Luer adapters.

OVERDOSE
 Not applicable if properly used in relation to intended use and in relation to catheter volume.

CONTRAINDICATIONS AND SIDE EFFECTS
 No contraindications. There are no known side effects if the solution is used as intended.

INTERACTIONS WITH OTHER SUBSTANCES OR MEDICAL PRODUCTS
 None when used for its intended purpose, i.e. for infusing bladder catheter balloons.

INSTRUCTIONS FOR USE
 Single-use device intended for use for a single patient.
 In order to ensure sterile application of the product, the healthcare professional must comply with the procedures of their institution, adopting an aseptic technique for the entire procedure.
 Check the validity of the expiry date indicated on the sealing side of the pouch.
 Open the pouch and take the syringe out. Check that the syringe cap is positioned properly and the solution is transparent. If necessary, place the syringe on a sterile field.
 Remove the closure cap, ensuring that no contact contamination of the syringe Luer cone occurs.
 Connect the syringe to the balloon inflation channel on the catheter, taking care to avoid any contamination through contact with the syringe Luer cone.
 Ensure that the quantity of solution used is equal to the volume of the balloon which it is to be inserted. Refer to the catheter manufacturer's instructions for correct identification of the catheter volume.
 Once the syringe has been connected to the catheter, slowly introduce into the catheter the amount of solution required to properly and effectively secure the balloon inside the bladder catheter, in accordance with the procedure adopted by your institution.
 Make sure that the volume of solution introduced is equal to the capacity of the balloon being filled, in order to avoid overfilling the balloon, as could result in bursting.
 After use, discard the empty syringe and any unused volumes of solution in accordance with national regulations.

SHELF LIFE
 3 years in unopened packaging.

STORAGE
 Store the syringes lightly closed, away from direct light and heat sources, at a temperature between +5 °C and +30 °C. Do not freeze. The expiry dates refers to the product in unopened packaging, correctly stored. Do not use the Medical Device after this date.

REPORTING OF SERIOUS INCIDENTS
 The user and/or patient is advised to report any serious incidents occurring in relation to the Medical Device to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

METHOD OF DISPOSAL
 The Medical Device is for single use. After use, dispose of the product in the appropriate hospital containers and according to the protocols of the health facility where the product is used or according to the regulations of your municipality. Do not release to the environment.

SYMBOLS USED ON THE LABEL	
	Do not reuse - Disposable
	Manufacturer
	Product code
	Production batch
	Expiry date
	Do not use if the pouch is open or damaged
	Consult instructions for use
	Latex-free
	Temperature limits for storage
	Protect from moisture
	Medical Device compliant with EU Regulation 2017/745

D.B.M. S.r.l. - Single member limited liability company
 Via dell'Industria, 2 - 23030 Tovo di Sant'Agata (SO) - Italy
 T +39 0342 771613 - info.DBM@numantec.eu - www.dmbitalia.it

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The volume reported in the label refers to the syringe filling volume. The graduated scale on the label of the immediate packaging (syringe) has no measuring function.

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